



More information on the website  
[mirror.radwag.com/en/info,w1,12Q](https://mirror.radwag.com/en/info,w1,12Q)

# Suspended Self-Centring Weighing Pan for WAY 500.5Y.KO, WAY 1.5Y.KO, WAY 2.5Y.KO Manual Mass Comparator

WX-001-0426



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

The suspended self-centring weighing pan is designed specifically for the WAY 5Y.KO manual mass comparators.

Installed It is installed during the device's production stage in place of the standard weighing pan, so it should be ordered when ordering a comparator (it is not sold as a separate product). The weighing pan has the following functions: placing several weights, which are the resultant of the reference standard, off-centre, positioning the centre of gravity of the test and reference standards always in the same place.



## Compatible with (Additional Fee)

WAY 5Y.KO Manual Mass Comparator