



More information on the website
mirror.radwag.com/en/info,w1,6NX

Cable PT0453

WX-013-0238



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

Compatible with (Additional Fee)

XA 5Y.M.A Microbalance
 XA 5Y.A Analytical Balance
 AS 5Y Analytical Balance
 XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance
 Automatic feeders
 XA 5Y.M Microbalance
 Receipt Printer