



More information on the website  
[mirror.radwag.com/en/info,w1,6NX](https://mirror.radwag.com/en/info,w1,6NX)

## PT0453 Cable

WX-013-0238



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

### Compatible with (Additional Fee)

XA 5Y.M.A Microbalance  
XA 5Y.A Analytical Balance  
AS 5Y Analytical Balance  
XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance  
Automatic feeders  
XA 5Y.M Microbalance  
Receipt Printer