



More information on the website  
[mirror.radwag.com/en/info,w1,6NX](https://mirror.radwag.com/en/info,w1,6NX)

## PT0453 Cable



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

### Compatible with

XA 5Y.M.A Microbalance  
XA 5Y.A Analytical Balance  
AS 5Y Analytical Balance  
XA 5Y.F Analytical Balance for Filters

XA 5Y Analytical Balance  
Automatic feeders  
XA 5Y.M Microbalance  
Receipt Printer