



More information on the website
mirror.radwag.com/us/info,w1,6NX

PT0453 Cable

WX-013-0238



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

Compatible with (Additional Fee)

XA 5Y.M.A Microbalance
XA 5Y.A Analytical Balance
AS 5Y Analytical Balance
Analytical Balance for Filters XA 5Y.F

XA 5Y Analytical Balance
Automatic feeders
XA 5Y.M Microbalance
Receipt Printer