




# RES Radwag Extended SI Module (5Y, CY10)

WX-016-0264

More information on the website  
[mirror.radwag.com/br/info,w1,8RT](http://mirror.radwag.com/br/info,w1,8RT)



The drawings, photos and graphics used are for illustrative purposes only.

## Visão geral

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi\* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

\* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

## Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P  
Microbalança para calibração de micropipetas MYA 5Y.P  
Analizador de umidade MA 5Y  
Balança de precisão 5Y.PM  
Microbalança para filtros MYA 5Y.F  
Microbalança XA 5Y.M.A

Comparador de massa manual APP 5Y.KO  
Medição de densidade  
Microbalança MYA 5Y  
Balança analítica XA 5Y  
Susceptômetros manuais - Magnetização  
Comparadores de massa manuais 5Y.PM.KB

!Terminal de pesagem PUE CY10  
Balanças automáticas para calibração de micropipetas  
Balança analítica XA 5Y.A  
Comparador de massa manual WAY 5Y.KO  
Comparadores de massa robóticos  
Ultra-Microbalança UYA 5Y  
Analizador de umidade MA 5Y.IC  
Balança analítica AS 5Y  
Comparador de massa automático AK-4  
Comparador de massa automático AKM-2  
Balança de precisão CY10  
Balança analítica para filtros XA 5Y.F  
Comparadores de massa automáticos UMA

Comparador de massa HRP 5Y.KO  
Balança de precisão PS 5Y  
Comparadores de massa automáticos de vácuo  
Analizador de umidade PMV 5Y  
Susceptômetros automáticos - Magnetização  
Comparador de massa manual XA 5Y.KB  
Comparador de massa HRP 5Y.KB  
Microbalança XA 5Y.M  
Comparador de massa manual UYA 5Y.KO  
Ultra-Microbalança para filtros UYA 5Y. F  
Nano-Comparador NANO.AK-4.500.5Y  
Microbalança para pesagem de stents XA 5Y.M.A.S