



More information on the website  
[mirror.radwag.com/it/info,w1,AFX](https://mirror.radwag.com/it/info,w1,AFX)

## Filter tray 160

WX-001-0405



The drawings, photos and graphics used are for illustrative purposes only.

### Caratteristiche generali

The tray collects particles that may fall off the filter. It is an additional protection of the weighing mechanism. The tray should be placed inside the filter chamber under the weighing pan on which the filter is weighed. It can be easily removed and cleaned.



### Compatibile con (Additional Fee)

Microbilancia per filtri MYA 5Y.F