

KIT 195 Density Determination Kit

WX-002-0021



More information on the website mirror.radwag.com/us/info,w1,FS4



The drawings, photos and graphics used are for illustrative purposes only.

Overview

Density determination kit for determining density of solids and liquids is designed for cooperation with laboratory balances featuring 195×195 mmweighing pan. This stainless steel kit is applicable for determining density of: solids; liquids. For testing granulates of density value higher or lower than 1 g/cm3 use KIT 85.

Compatible with (Additional Fee)

WLC X7 Precision balance WLC Precision Balance

WLC C/2 Precision balance PS R2 Precision Balance

Device dimensions W x D x H

(3,9)


12 ±


.....


23±0,3

 \triangle


105,6







KIT 128



Granulte density determination KIT

KIT 195.M

R125

R12,5

00

KIT 85

