



More information on the website
mirror.radwag.com/us/info,w1,GWF

Automatic Chamber for filter weighing

WX-002-0135



The drawings, photos and graphics used are for illustrative purposes only.

Overview

The chamber for filter weighing is characterized by high ingress protection. Large open-work weighing pan enables precise weighing of filters of various dimensions (max $\varnothing 100$). The chamber is made of stainless steel and is resistant to chemicals, colour change and scratches. The new system of chamber disassembly facilitates maintenance activities. The top lid of the chamber has been equipped with a glass opening allowing the operator to check whether the filter remains in its start position after chamber closing.

The chamber for filter weighing may be used in MYA 5Y microbalance and UYA 5Y ultra-microbalance after disassembly of a regular weighing chamber.

When placing an order for a chamber for filter weighing, a serial number of the scale is required.

Compatible with (Additional Fee)

MYA 5Y.P Microbalance for pipettes calibration
UYA 5Y Ultra-Microbalance

MYA 5Y Microbalance