



More information on the website
mirror.radwag.com/en/info,w1,IRH

20 μ l ÷ 200 μ l - RW8-105-20-9 Automatic Pipettes

PA-203-0005



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters

| | |
|------------------|------------------|
| Volume | 20 – 200 μ l |
| Systematic error | 2.0% / 0.6% |
| Random error | 0.8% / 0.15% |
| Tips | 300 μ l |



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration

Automatic Balances for Pipettes Calibration