



More information on the website
mirror.radwag.com/br/info,w1,OC3

Cable PT0455

WX-013-0239



The drawings, photos and graphics used are for illustrative purposes only.

Visão geral

PT0455 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a RADWAG RTP printer. Its standard length is 1,8m.

Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P
 Microbalança para calibração de micropipetas MYA 5Y.P
 Balança de precisão 5Y.PM
 Microbalança para filtros MYA 5Y.F
 Microbalança XA 5Y.M.A
 Balanças automáticas para calibração de micropipetas
 Comparador de massa manual WAY 5Y.KO
 Comparadores de massa robóticos
 Ultra-Microbalança UYA 5Y
 Comparador de massa automático AK-4
 Comparador de massa automático AKM-2
 Balança analítica para filtros XA 5Y.F
 Comparadores de massa automáticos UMA
 Comparador de massa manua APP 5Y.KO
 Medição de densidade

Microbalança MYA 5Y
 Susceptômetros manuais - Magnetização
 Comparadores de massa manuais 5Y.PM.KB
 Comparador de massa HRP 5Y.KO
 Comparadores de massa automáticos de vácuo
 Analisador de umidade PMV 5Y
 Susceptômetros automáticos - Magnetização
 Comparador de massa manual XA 5Y.KO
 Comparador de massa HRP 5Y.KB
 Microbalança XA 5Y.M
 Comparador de massa manual UYA 5Y.KO
 Ultra-Microbalança para filtros UYA 5Y. F
 Nano-Comparador NANO.AK-4.500.5Y
 Microbalança para pesagem de stents XA 5Y.M.A.S
 Impressora de recibos