



More information on the website  
[mirror.radwag.com/br/info,w1,OC3](http://mirror.radwag.com/br/info,w1,OC3)

## Cable PT0455

WX-013-0239



The drawings, photos and graphics used are for illustrative purposes only.

## Visão geral

PT0455 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a RADWAG RTP printer. Its standard length is 1,8m.

## Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P  
 Microbalança para calibração de micropipetas MYA 5Y.P  
 Balança de precisão 5Y.PM  
 Microbalança para filtros MYA 5Y.F  
 Microbalança XA 5Y.M.A  
 Balanças automáticas para calibração de micropipetas  
 Comparador de massa manual WAY 5Y.KO  
 Comparadores de massa robóticos  
 Ultra-Microbalança UYA 5Y  
 Comparador de massa automático AK-4  
 Comparador de massa automático AKM-2  
 Balança analítica para filtros XA 5Y.F  
 Comparadores de massa automáticos UMA  
 Comparador de massa manua APP 5Y.KO

Microbalança MYA 5Y  
 Susceptômetros manuais - Magnetização  
 Comparadores de massa manuais 5Y.PM.KB  
 Comparador de massa HRP 5Y.KO  
 Comparadores de massa automáticos de vácuo  
 Analisador de umidade PMV 5Y  
 Susceptômetros automáticos - Magnetização  
 Comparador de massa manual XA 5Y.KO  
 Comparador de massa HRP 5Y.KB  
 Microbalança XA 5Y.M  
 Comparador de massa manual UYA 5Y.KO  
 Ultra-Microbalança para filtros UYA 5Y. F  
 Nano-Comparador NANO.AK-4.500.5Y  
 Microbalança para pesagem de stents XA 5Y.M.A.S

