



More information on the website
mirror.radwag.com/en/info,w1,OC3

Cable PT0455

WX-013-0239



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0455 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a RADWAG RTP printer. Its standard length is 1,8m.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration
UMA 5Y.F Automatic Weighing System
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
Automatic Balances for Pipettes Calibration
WAY 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
UYA 5Y Ultra-Microbalance
AK-4 Automatic Mass Comparator
AKM-2 Automatic Mass Comparator
XA 5Y.F Analytical Balance for Filters
UMA Automatic Mass Comparator

Manual susceptometers - Magnetization
5Y.PM.KB Manual Mass Comparator
HRP 5Y.KO Mass Comparator
Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
Automatic susceptometers - Magnetization
XA 5Y.KO Manual Mass Comparator
HRP 5Y.KB Mass Comparator
RMC F Robotic Weighing System
RB F Robotic Weighing System
XA 5Y.M Microbalance
UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
Nano-Comparator NANO.AK-4.500.5Y

APP 5Y.KO Manual Mass Comparator
Density Measurement
MYA 5Y Microbalance

XA 5Y.M.A.S Microbalance for stents weighing
Receipt Printer