



More information on the website  
[mirror.radwag.com/en/info,w1,OC3](http://mirror.radwag.com/en/info,w1,OC3)

# Cable PT0455

WX-013-0239



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

PT0455 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a RADWAG RTP printer. Its standard length is 1,8m.

## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 UMA 5Y.F Automatic Weighing System  
 5Y.PM Precision Balance  
 MYA 5Y.F Microbalance for Filters  
 XA 5Y.M.A Microbalance  
 Automatic Balances for Pipettes Calibration  
 WAY 5Y.KO Manual Mass Comparator  
 Robotic Mass Comparators  
 UYA 5Y Ultra-Microbalance  
 AK-4 Automatic Mass Comparator  
 AKM-2 Automatic Mass Comparator  
 XA 5Y.F Analytical Balance for Filters  
 UMA Automatic Mass Comparator

Manual susceptometers - Magnetization  
 5Y.PM.KB Manual Mass Comparator  
 HRP 5Y.KO Mass Comparator  
 Automatic Vacuum Mass Comparators  
 PMV 5Y Moisture Analyzer  
 Automatic susceptometers - Magnetization  
 XA 5Y.KO Manual Mass Comparator  
 HRP 5Y.KB Mass Comparator  
 RMC F Robotic Weighing System  
 RB F Robotic Weighing System  
 XA 5Y.M Microbalance  
 UYA 5Y.KO Manual Mass Comparator  
 UYA 5Y.F Ultra-Microbalance for Filters  
 Nano-Comparator NANO.AK-4.500.5Y

APP 5Y.KO Manual Mass Comparator  
Density Measurement  
MYA 5Y Microbalance

XA 5Y.M.A.S Microbalance for stents weighing  
Receipt Printer