



More information on the website
mirror.radwag.com/us/info,w1,QLH

5 μ l ÷ 50 μ l - RW8-103-20-9 Automatic Pipettes

PA-203-0003



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| | |
|------------------|----------------|
| Volume | 5 – 50 μ l |
| Systematic error | 2,0% / 0,6% |
| Random error | 2,0% / 0,3% |
| Tips | 300 μ l |



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration

Automatic Balances for. Pipettes Calibration
Pipette Tips