



More information on the website  
[mirror.radwag.com/it/info,w1,QNI](https://mirror.radwag.com/it/info,w1,QNI)

## Voluminous sample container



The drawings, photos and graphics used are for illustrative purposes only.

### Caratteristiche generali

The set is an ideal solution for drying large volume samples, e.g. textiles, or products that placed on the weighing pan tend to increase their volume. The container allows to enclose the sample within a limited space wherein the possibility to determine the moisture content value remains unaffected.

In the case of large-volume samples, the kit prevents the sample from coming into contact with solid elements of the moisture analyzer (such as the heating element or its cover) and protects it from damage. Closing a large-volume, low-density sample in the container also reduces the risk of errors resulting from the effect of the sample being raised by warm air.

The accessory gives the user full control of the sample and an error-free result.

The set consists of 2 components made of steel. It is compatible with all moisture analyzers.

### Compatibile con

Termobilancia MA X7.IC.A  
Termobilancia MA 5Y  
Termobilancia MA X2.IC.A

Termobilancia MA X7.A  
Termobilancia MA 5Y.IC  
Termobilancia MA R

