




Fingerprint Reader

WX-006-0048

More information on the website
mirror.radwag.com/en/info,w1,Q01



The drawings, photos and graphics used are for illustrative purposes only.

Overview

An accessory for the ELLIPSIS laboratory balances, designed to read fingerprints in order to log into the balance. Simply touch the reader to log in. The fingerprint reader is connected to the ELLIPSIS balance via the USB port.



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MYA 5Y.P Microbalance for pipettes calibration
AK-5/100.STENT Automatic Balance for Stents
UMA 5Y.F Automatic Weighing System
MA 5Y Moisture Analyzer
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
CY10 Scales for PGC

APP 5Y.KO Manual Mass Comparator
Robotic Mass Comparators
Density Measurement
MYA 5Y Microbalance
XA 5Y Analytical Balance
Manual susceptometers - Magnetization
5Y.PM.KB Manual Mass Comparator
HRP 5Y.KO Mass Comparator
PS 5Y Precision Balance

PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
Automatic Balances for. Pipettes Calibration
XA 5Y.A Analytical Balance
WAY 5Y.KO Manual Mass Comparator
UYA 5Y Ultra-Microbalance
MA 5Y.IC Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator
PS 5Y Balances for PGC
AKM-2 Automatic Mass Comparator
CY10 Precision Balance
XA 5Y.F Analytical Balance for Filters
UMA Automatic Mass Comparator

Automatic Vacuum Mass Comparators
PMV 5Y Moisture Analyzer
Automatic susceptometers - Magnetization
XA 5Y.KO Manual Mass Comparator
HRP 5Y.KB Mass Comparator
XA 5Y.M Microbalance
RMC F Robotic Weighing System
RB F Robotic Weighing System
UYA 5Y.KO Manual Mass Comparator
UYA 5Y.F Ultra-Microbalance for Filters
Nano-Comparator NANO.AK-4.500.5Y
XA 5Y.M.A.S Microbalance for stents weighing
AK-6 5Y.F Automatic Weighing System