




# Protocolo de comunicación Radwag (X2)



The drawings, photos and graphics used are for illustrative purposes only.

## Descripción

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi\* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

\* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

## Trabaja con

Analizador de Humedad  
MA X7.IC.A  
Analizadores de HumedadMA X2.IC.A  
Balanza analítica AS X7  
Balanza de precisión PS X2  
Balanzas analíticas AS X2 PLUS

Balanzas PS X7 para CCE  
Balanzas PS X2 para CCE  
Analizadores de HumedadMA X2.A  
Analizador de Humedad  
MA X7.A  
Balanza de precisión PS X7