



More information on the website
mirror.radwag.com/it/info,w1,W13

Protocollo di comunicazione Radwag (X2)



The drawings, photos and graphics used are for illustrative purposes only.

Caratteristiche generali

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

Compatibile con

Termobilancia MA X7.IC.A
Termobilance MA X2.IC.A
AS X7 Bilancia analitica
Bilance di precisione PS X2
Bilance analitiche AS X2 PLUS

Bilance PS X7 per controllo CPC
Bilance PS X2 per controllo CPC
Termobilance MA X2.A
Termobilancia MA X7.A
Bilancia di precisione PS X7