




More information on the website  
[mirror.radwag.com/it/info,w1,W13](http://mirror.radwag.com/it/info,w1,W13)

# Protocollo di comunicazione Radwag (X2)

WX-016-0265



The drawings, photos and graphics used are for illustrative purposes only.

## Caratteristiche generali

Character communication protocol is designed for communication between RADWAG balance and external device via RS232, Ethernet or Wi-Fi\* interfaces. The protocol consists of commands sent from the external device to the balance and responses sent from the balance to the device. Responses are sent from the balance each time a command is received as a response to the command. The commands that make up the communication protocol allow you to obtain information about the status of the balance and ways to influence its operation, for example, they allow you to receive weighing results from the balance, control the display, etc.

\* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.

## Compatibile con (Additional Fee)

Termobilancia MA X7.IC.A  
AS X7 Bilancia analitica  
Bilance PS X7 per controllo CPC

Termobilancia MA X7.A  
Bilancia di precisione PS X7