



More information on the website  
[mirror.radwag.com/en/info,w1,01Z](https://mirror.radwag.com/en/info,w1,01Z)

# Knob weights, Set (1 mg - 10 kg), Plastic Box

OK-502-0021



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 21 111.11 g     |
| Nominal value           | 1 mg – 10 kg    |
| Quantity                | 29 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes