



More information on the website
mirror.radwag.com/en/info,w1,07X

IM02.1 (MediaBox) Communication Module

WX-016-0303



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Protection class	IP 43
Communication interface	
Communication interface	3xRS232, 4IN/4OUT, Virtual Com
Environmental conditions	
Relative humidity	10% – 85% RH no condensation
Physical parameters	
Device dimensions W x D x H	100×100×37 mm



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration

CY10 Precision Balance

MYA 5Y.P Microbalance for pipettes calibration
 5Y.PM Precision Balance
 MYA 5Y.F Microbalance for Filters
 XA 5Y.M.A Microbalance
 CY10 Scales for PGC
 PUE CY10 Weighing Terminal
 CY10 Multifunctional Scale
 XA 5Y.A Analytical Balance
 UYA 5Y Ultra-Microbalance

XA 5Y.F Analytical Balance for Filters
 MYA 5Y Microbalance
 XA 5Y Analytical Balance
 PMV 5Y Moisture Analyzer
 XA 5Y.M Microbalance
 RS 232, RS 485 cables
 UYA 5Y.F Ultra-Microbalance for Filters
 Additional modules
 XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions W x D x H

