



More information on the website  
[mirror.radwag.com/en/info,w1,07X](http://mirror.radwag.com/en/info,w1,07X)

# IM02.1 (MediaBox) Communication Module

WX-016-0303



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Protection class	IP 43
Communication interface	
Communication interface	3xRS232, 4IN/4OUT, Virtual Com
Environmental conditions	
Relative humidity	10% – 85% RH no condensation
Physical parameters	
Device dimensions W x D x H	100×100×37 mm



## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration

CY10 Precision Balance

MYA 5Y.P Microbalance for pipettes calibration  
 5Y.PM Precision Balance  
 MYA 5Y.F Microbalance for Filters  
 XA 5Y.M.A Microbalance  
 CY10 Scales for PGC  
 PUE CY10 Weighing Terminal  
 CY10 Multifunctional Scale  
 XA 5Y.A Analytical Balance  
 UYA 5Y Ultra-Microbalance

XA 5Y.F Analytical Balance for Filters  
 MYA 5Y Microbalance  
 XA 5Y Analytical Balance  
 PMV 5Y Moisture Analyzer  
 XA 5Y.M Microbalance  
 RS 232, RS 485 cables  
 UYA 5Y.F Ultra-Microbalance for Filters  
 Additional modules  
 XA 5Y.M.A.S Microbalance for stents weighing

## Device dimensions W x D x H

