



More information on the website
mirror.radwag.com/us/info,w1,07X

IM02-1 (MediaBox) Communication Module



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Communication interface	RS232, 4IN/4OUT, Virtual Com
Physical parameters	
Device dimensions	100×100×30,5 mm




Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 5Y.PM Precision balance
 MYA 5Y.F Microbalance for Filters
 XA 5Y.M.A Microbalance
 CY10 Scales for PGC
 PUE CY10 Weighing Terminal
 CY10 Multifunctional Scale
 XA 5Y.A Analytical Balance
 UYA 5Y Ultra-Microbalance

CY10 Precision balance
 Analytical Balance for Filters XA 5Y.F
 MYA 5Y Microbalance
 XA 5Y Analytical Balance
 PMV 5Y Moisture Analyzer
 XA 5Y.M Microbalance
 RS 232, RS 485 cables
 UYA 5Y.F Ultra-Microbalance for Filters
 XA 5Y.M.A.S Microbalance for stents weighing

Device dimensions



100×100×30,5 mm