



More information on the website
mirror.radwag.com/en/info,w1,0F6

Transport Case for PS, WLC ((X7, A1 and A2) Balances

WX-002-0049



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Physical parameters

Packaging dimensions W x D x H

520×370×145 mm



Compatible with (Additional Fee)

PS X7 Balances for PGC
PS R2 Precision Balance
WLC X7 Precision Balance
WLC Precision Balance

WLC C/2 Precision Balance
PS R1 Precision Balance
PS X7 Precision Balance

Device dimensions W x D x H

