



More information on the website  
[mirror.radwag.com/en/info,w1,0F6](https://mirror.radwag.com/en/info,w1,0F6)

# Transport Case for PS, WLC ((X7, A1 and A2) Balances

WX-002-0049



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Physical parameters

Packaging dimensions W x D x H

520×370×145 mm



## Compatible with (Additional Fee)

PS X7 Balances for PGC  
PS R2 Precision Balance  
WLC X7 Precision Balance  
WLC Precision Balance

PS X7 Precision Balance  
WLC C/2 Precision Balance  
PS R1 Precision Balance

## Device dimensions W x D x H

