



More information on the website  
[mirror.radwag.com/en/info,w1,0YZ](https://mirror.radwag.com/en/info,w1,0YZ)

## BT 10 µl Tips - 1000 pcs

PT-110-0001



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

OMNITIP™ BT 10 µl tips for pipettes are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

**Capacity: 10 µl Colour: transparent Set contains 1000 tips in a plastic bag with a zipper closure**

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

### Compatible with (Additional Fee)

