



More information on the website  
[mirror.radwag.com/us/info,w1,0YZ](https://mirror.radwag.com/us/info,w1,0YZ)

## BT 10 $\mu$ l Tips - 1000 pcs

PT-110-0001



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

OMNITIP™ - Bulk BT 10, transparent. 10  $\mu$ l, 1000 tips in bulk.

### Compatible with (Additional Fee)

Workstation for pipettes calibration

Automatic Variable-Volume Pipettes