




More information on the website
mirror.radwag.com/en/info,w1,1KQ

RAD Key

WX-010-0005



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Components and software

Operating system	Windows 7, 8, 10
------------------	------------------

Compatible with (Additional Fee)

H315 HR Waterproof Platform Scale
 MA X7.IC.A Moisture Analyzer
 HX7 H Stainless steel Multifunctional Scale
 HX7.4N H Multifunctional Stainless Steel Ramp Scale
 HX7.4C Multifunctional Platform Scale
 HX7 HR Multifunctional Scale
 XA 5Y.M.A.P Microbalance for pipettes calibration
 WPY Multifunctional Scale
 MYA 5Y.P Microbalance for pipettes calibration
 C32 Multifunctional Scale
 Stainless Steel HRP EX.H Platforms
 HX5.EX-1.4P2 H Stainless Steel Beam Scales
 HY10.HRP.H Stainless steel High Resolution Scale
 MA 5Y Moisture Analyzer
 C315 Platform scale
 MYA 5Y.F Microbalance for Filters
 HX5.EX-1 HR Stainless Steel 1-Load-Cell Platform Scales
 HX7.4P2 H Multifunctional Stainless Steel Beam Scales

XA 5Y.A Analytical Balance
 HX7.4 H Multifunctional Stainless Steel Platform Scale
 C315.4P pallet scales
 HX7.4P2 Multifunctional Beam Scales
 WAY 5Y.KO Manual Mass Comparator
 UYA 5Y Ultra-Microbalance
 MA 5Y.IC Moisture Analyzer
 HX5.EX-1 C 1-Load-Cell Platform Scales
 MA R Moisture Analyzer
 HX5.EX-1 H Waterproof Platform Scales
 AS 5Y Analytical Balance
 AK-4 Automatic Mass Comparator
 PS 5Y Precision Balance
 PS 5Y Balances for PGC
 AKM-2 Automatic Mass Comparator
 PS R2.H Precision Balance
 PUE HX5.EX Weighing Terminal
 HY10.H Scales for PGC

TMX H Stainless steel Multifunctional Scale	HX5.EX-1.4 H Stainless Steel 4-Load-Cell Platform Scales
WLY Scales for PGC	WTC Precision Balance
AS X7 Analytical Balance	CY10 Precision Balance
HRP Platforms High Capacity	H315 4C Platform Scale
PS X7 Balances for PGC	C315.OW Personal Scale
TMX Multifunctional Scale	XA 5Y.F Analytical Balance for Filters
C315.D Baby Scale	WLY Precision Balance
MWMH Series Weighing Modules	UMA Automatic Mass Comparator
XA 5Y.M.A Microbalance	PS X7 Precision Balance
WLC Precision Balance	APP 5Y.KO Manual Mass Comparator
PS R2 Precision Balance	H315 4P.H Stainless Steel Pallet Scales
C315.4P2 Beam Scales	Robotic Mass Comparators
H315 4H Stainless Steel Platform Scale	MW-01 Mass Converters
H315 4P2.H Stainless steel Beam Scales	WLC C/2 Precision Balance
WLC X7 Precision Balance	MYA 5Y Microbalance
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale	PS R1 Precision Balance
HX5.EX-1.4P2 C Beam Scales	XA 5Y Analytical Balance
CY10 Scales for PGC	PUE HX7 Weighing Terminal
Stainless steel HRP.H platforms	5Y.PM.KB Manual Mass Comparator
H315 H Waterproof Platform Scale	AS R2 PLUS Analytical Balance
PUE CY10 Weighing Terminal	PUE C32 Weighing Terminal
CY10 Multifunctional Scale	PUE 7.1 Weighing Terminal
HX7.4P H Multifunctional Stainless Steel Pallet Scales	PUE 5 Weighing Terminal
HX5.EX-1.4N Stainless Steel Ramp Scales	HRP 5Y.KO Mass Comparator
HY10 Multifunctional Scale	Postal Scales for Packages
HX5.EX-1.4 C 4-Load-Cell Platform Scales	PMV 5Y Moisture Analyzer
HY10 H Stainless Steel Multifunctional Scale	HX5.EX-1.4P C Pallet Scales
H315.4 H/Z Stainless Steel Platform Scale, pit version	XA 5Y.KO Manual Mass Comparator
H315 4N.H Stainless Steel Ramp Scale	HRP 5Y.KB Mass Comparator
HRP EX Platforms High Capacity	XA 5Y.M Microbalance
C315.4 Platform Scale	PUE C315 Weighing Terminal
HX5.EX-1.4P H Stainless Steel Pallet Scales	PUE HY10 Weighing Terminal
HY10.HRP High Resolution Scale	MW-04 Mass Converters
C32.PM Precision Balance	UYA 5Y.KO Manual Mass Comparator
PUE 7.1.HRP High Resolution Scale	AS R1 PLUS Analytical Balance
Automatic Balances for Pipettes Calibration	UYA 5Y.F Ultra-Microbalance for Filters
HX7 Multifunctional Scale	C315.4B Bed Scale
MWLH Series Weighing Modules	PUE H315 Weighing Terminal
HX5.EX-1.4 H/Z Stainless Steel Platform Scales, Pit Version	HX7.4P Multifunctional pallet scales
MA X7.A Moisture Analyzer	Nano-Comparator NANO.AK-4.500.5Y
HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version	XA 5Y.M.A.S Microbalance for stents weighing
C315.O Personal Scale	MAS.1.Y Modules
H315 4N.H.LD Stainless Steel Ramp Scale	C315.8B Ramps Bed Scale
WPY Scales for PGC	