



More information on the website  
[mirror.radwag.com/en/info,w1,1N2](https://mirror.radwag.com/en/info,w1,1N2)

# Cylinder Weights, Set (1 g - 200 g), Plastic Box

OK-102-0049



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	610 g
Nominal value	1 – 200 g
Quantity	11 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes