



More information on the website  
[mirror.radwag.com/tr/info,w1,10F](http://mirror.radwag.com/tr/info,w1,10F)

## UL-04-1 (AS) Tüp Tutucu

WX-001-0033



The drawings, photos and graphics used are for illustrative purposes only.

### Teknik Veriler

#### Fiziksel parametreler

Ambalaj boyutları G x D x Y	110×102×56 mm
Net ağırlık	103 g
Brüt ağırlık	145 g

### İle uyumlu (Additional Fee)

AS X7 Analitik Teraziler  
AS 5Y Analitik Teraziler

AS R2 PLUS Analitik Teraziler  
AS R1 PLUS Analitik Teraziler

## Cihaz boyutları G x D x Y

