



More information on the website  
[mirror.radwag.com/de/info,w1,23D](https://mirror.radwag.com/de/info,w1,23D)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 g - 1 kg), Holz-Etui

OK-401-0075



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Messtechnische Parameter

OIML-Klasse	M1
Gewicht	2 110 g
Nominalwert	1 g - 1 kg
Anzahl	13 St.

## Kompatibel mit (Additional Fee)

Pinsel