



More information on the website  
[mirror.radwag.com/en/info,w1,23D](https://mirror.radwag.com/en/info,w1,23D)

# M1 Mass Standard, Knob Weights, Set (1 g - 1 kg), Wooden Box

OK-401-0075



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	M1
Weight	2 110 g
Nominal value	1 g – 1 kg
Quantity	13 pcs

## Compatible with (Additional Fee)

Brushes