



More information on the website
mirror.radwag.com/it/info,w1,23D

Pesi standard M1, cilindrici, set (1 g - 1 kg), cassa di legno

OK-401-0075



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici

| | |
|-----------------|------------|
| Classe OIML | M1 |
| Peso | 2 110 g |
| Valore nominale | 1 g - 1 kg |
| Numero | 13 pcs |

Compatibile con (Additional Fee)

Pennelli