



More information on the website  
[mirror.radwag.com/de/info,w1,2BV](https://mirror.radwag.com/de/info,w1,2BV)

# Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 g - 200 g), Holz-Etui

OK-401-0073



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

### Messtechnische Parameter

|             |           |
|-------------|-----------|
| OIML-Klasse | M1        |
| Gewicht     | 610 g     |
| Nominalwert | 1 – 200 g |
| Anzahl      | 11 St.    |

## Kompatibel mit (Additional Fee)

Pinsel