



More information on the website
mirror.radwag.com/en/info,w1,2C8

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 2 kg), Wooden Box

OK-301-0042



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	6 111.11 g
Nominal value	1 mg – 2 kg
Quantity	27 pcs

Compatible with (Additional Fee)

Tweezers

Brushes