



More information on the website
mirror.radwag.com/en/info,w1,2FV

Knob weights, Set (1 g - 50 g), Wooden Box

OK-501-0027



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | E1 |
| Construction | |
| Material | Stainless steel |
| Weight | 110 g |
| Nominal value | 1 – 50 g |
| Quantity | 8 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes