



More information on the website
mirror.radwag.com/de/info,w1,2XC

Prüfgewichte der E1 Klasse - Knopfgewichte, Satz (1 mg - 100 g), Kunststoff-Etui

OK-502-0015



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Gewicht	211,11 g
Nominalwert	1 mg – 100 g
Anzahl	21 St.

Kompatibel mit (Additional Fee)

Pinzetten

Pinzel