



More information on the website  
[mirror.radwag.com/en/info,w1,2XC](http://mirror.radwag.com/en/info,w1,2XC)

# Knob weights, Set (1 mg - 100 g), Plastic Box

OK-502-0015



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	211.11 g
Nominal value	1 mg – 100 g
Quantity	21 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes