



More information on the website  
[mirror.radwag.com/en/info,w1,3BF](https://mirror.radwag.com/en/info,w1,3BF)

# Stainless Steel Wire, Set (1 mg - 500 mg), Plastic Box

OK-502-0002



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	1.11 g
Nominal value	1 mg - 500 mg
Weights quantity	12

## Compatible with (Additional Fee)

Tweezers

Brushes