



More information on the website  
[mirror.radwag.com/en/info,w1,3DP](https://mirror.radwag.com/en/info,w1,3DP)

# F1 Mass Standard - cylindrical weight - 20

**g**

OP-203-0116



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Nominal value           | 20 g            |

## Compatible with (Additional Fee)

Tweezers

Brushes