



More information on the website
mirror.radwag.com/us/info,w1,3KF

ZR-04/1 Power adapter

WX-015-0329



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Power supply	Type-C: 5V/3A, 9V/3A, 12V/3A, 15V/3A, 20V/3A 60W (max) (ZR-04)
Physical parameters	
Device dimensions W x D x H	153×69×42 mm
Packaging dimensions W x D x H	230×190×80 mm
Net weight	645 g
Gross weight	900 g



Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
 MYA 5Y.P Microbalance for pipettes calibration
 C32 Multifunctional Scale
 C315 Platform scale
 5Y.PM Precision balance
 MYA 5Y.F Microbalance for Filters

AS 5Y Analytical Balance
 PS 5Y Balances for PGC
 PS R2.H Precision Balance
 WTC Precision Balance
 CY10 Precision balance
 C315.OW Personal scale

WLY Scales for PGC
AS X7 Analytical Balance
PS X7 Balances for PGC
C315.D Baby scale
XA 5Y.M.A Microbalance
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
WLC X7 Precision balance
CY10 Scales for PGC
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
C315.4 Platform Scale
C32.PM Precision Balance
PS X7 Precision Balance
WPY Scales for PGC
XA 5Y.A Analytical Balance
C315.4P pallet scales
UYA 5Y Ultra-Microbalance

Analytical Balance for Filters XA 5Y.F
WLY Precision Balance
WLC C/2 Precision balance
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PS 5Y Precision Balance
XA 5Y.M Microbalance
PUE C315 Terminal
AS R1 PLUS Analytical Balance
UYA 5Y.F Ultra-Microbalance for Filters
C315.4B Bed scale
C315.K Chair scale
XA 5Y.M.A.S Microbalance for stents weighing
C315.8B Ramps bed scale