



More information on the website  
[mirror.radwag.com/us/info,w1,3KF](https://mirror.radwag.com/us/info,w1,3KF)

# ZR-04/1 Power adapter

WX-015-0329



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Electrical parameters	
Power supply	Type-C: 5V/3A, 9V/3A, 12V/3A, 15V/3A, 20V/3A 60W (max) (ZR-04)
Physical parameters	
Device dimensions W x D x H	153×69×42 mm
Packaging dimensions W x D x H	230×190×80 mm
Net weight	645 g
Gross weight	900 g



## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 C32 Multifunctional Scale  
 C315 Platform scale  
 5Y.PM Precision balance

PS 5Y Precision Balance  
 PS 5Y Balances for PGC  
 PS R2.H Precision Balance  
 WTC Precision Balance  
 CY10 Precision balance

MYA 5Y.F Microbalance for Filters  
WLY Scales for PGC  
AS X7 Analytical Balance  
PS X7 Balances for PGC  
C315.D Baby scale  
XA 5Y.M.A Microbalance  
WLC Precision Balance  
PS R2 Precision Balance  
C315.4P2 Beam Scales  
WLC X7 Precision balance  
CY10 Scales for PGC  
PUE CY10 Weighing Terminal  
CY10 Multifunctional Scale  
C315.4 Platform Scale  
C32.PM Precision Balance  
WPY Scales for PGC  
XA 5Y.A Analytical Balance  
C315.4P pallet scales  
UYA 5Y Ultra-Microbalance  
AS 5Y Analytical Balance

C315.OW Personal scale  
Analytical Balance for Filters XA 5Y.F  
WLY Precision Balance  
PS X7 Precision Balance  
WLC C/2 Precision balance  
MYA 5Y Microbalance  
PS R1 Precision Balance  
XA 5Y Analytical Balance  
AS R2 PLUS Analytical Balance  
PUE C32 Weighing Terminal  
PUE 7.1 Weighing Terminal  
XA 5Y.M Microbalance  
PUE C315 Terminal  
AS R1 PLUS Analytical Balance  
UYA 5Y.F Ultra-Microbalance for Filters  
C315.4B Bed scale  
C315.K Chair scale  
XA 5Y.M.A.S Microbalance for stents weighing  
C315.8B Ramps bed scale