



More information on the website  
[mirror.radwag.com/it/info,w1,400](http://mirror.radwag.com/it/info,w1,400)

# Pesi standard F1 - cilindrico - 100 g

OP-203-0104



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Classe OIML	F1
-------------	----

### Costruzione

Materiale	Stainless steel
-----------	-----------------

Valore nominale	100 g
-----------------	-------

## Compatibile con (Additional Fee)

Pinzette

Pennelli