



More information on the website  
[mirror.radwag.com/en/info,w1,4EU](https://mirror.radwag.com/en/info,w1,4EU)

# M1 Mass Standard, Knob Weights, Set (1 g - 5 kg), Wooden Box

OK-401-0077



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	11 110 g
Nominal value	1 g – 5 kg
Quantity	16 pcs

## Compatible with (Additional Fee)

Brushes