



More information on the website
mirror.radwag.com/it/info,w1,4EU

Pesi standard M1, cilindrici, set (1 g - 5 kg), cassa di legno

OK-401-0077



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri metrologici

| | |
|-----------------|------------|
| Classe OIML | M1 |
| Peso | 11 110 g |
| Valore nominale | 1 g - 5 kg |
| Numero | 16 pcs |

Compatibile con (Additional Fee)

Pennelli