



More information on the website  
[mirror.radwag.com/de/info,w1,57D](http://mirror.radwag.com/de/info,w1,57D)

# Prüfgewichte der E1 Klasse - Knopfgewichte, Satz (1 mg - 100 g), Holz- Etui

OK-501-0038



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Gewicht	211,11 g
Nominalwert	1 mg – 100 g
Anzahl	21 St.

## Kompatibel mit (Additional Fee)

Pinzetten

Pinsel