



More information on the website  
[mirror.radwag.com/en/info,w1,5AU](https://mirror.radwag.com/en/info,w1,5AU)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 5 kg), Wooden Box

OK-301-0033



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	11 110 g
Nominal value	1 g – 5 kg
Quantity	16 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes