



More information on the website  
[mirror.radwag.com/br/info,w1,6NX](https://mirror.radwag.com/br/info,w1,6NX)

## Cable PT0453

WX-013-0238



The drawings, photos and graphics used are for illustrative purposes only.

### Visão geral

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

### Compatível com (Additional Fee)

Microbalança XA 5Y.M.A  
Balança analítica XA 5Y.A  
Balança analítica AS 5Y  
Balança analítica para filtros XA 5Y.F

Balança analítica XA 5Y  
!Alimentadores automáticos  
Microbalança XA 5Y.M  
Impressora de recibos