

PT0453 Cable

WX-013-0238





The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0453 Cable is used to connect the 5Y series balance or the PUE CY10 terminal to a PA-04/H Automatic Feeder. Its standard length is 2 m.

Compatible with (Additional Fee)

XA 5Y.M.A Microbalance XA 5Y.A Analytical Balance AS 5Y Analytical Balance Analytical Balance for Filters XA 5Y.F XA 5Y Analytical Balance Automatic feeders XA 5Y.M Microbalance Receipt Printer